

Street Address:

AOCS, 2710 S. Boulder Drive Urbana, IL 61802-6996 USA

Phone: +1-217-359-2344 **Fax:** +1-217-351-8091 **E-Mail:** CRM@aocs.org **Web:** www.aocs.org

Certified Reference Materials AOCS 0610-A6

Report of the certification process for

FG72

Soybean Certified Reference Materials
Sixth Batch

Denise Williams Technical Services Manager Scott Bloomer Technical Service Director



Attention

BASF Agricultural Solutions Seed US LLC ("BASF") has produced reference material in support of the FG72 soybean event. BASF has transferred ownership of the FG72 soybean event to M.S. Technologies, L.L.C. ("MS Tech"), and pursuant to a service agreement between MS Tech and Syngenta Crop Protection, LLC, communications related to the reference material should be directed to Syngenta Crop Protection, 9 Davis Drive, Research Triangle Park, NC 27709-2252, USA

Legal Notice

Neither AOCS nor any person acting on behalf of AOCS is responsible for the use which might be made of the following information.

AOCS Mission Statement:

AOCS advances the science and technology of oils, fats, proteins, surfactants and related materials, enriching the lives of people everywhere.

More information regarding AOCS is available at http://www.aocs.org

Contents

Abstract	4
Acknowledgements	5
Glossary	6
Introduction	8
Materials and Methods	8
Stability	9
Results and Discussion	10
Sample Homogeneity and Prepared Sample Verification	10
References	12

Abstract

This report describes the preparation and certification of the soybean CRM AOCS 0610-A6 produced by AOCS Technical Services in 2021. The CRMs have been prepared according to ISO 17034:2016 and are intended to serve as control material for third party testing of soybean for transformation events. The presence of FG72 in the soybean was verified using event-specific, qualitative PCR analysis by FoodChain ID, Fairfield, IA (an ISO 17025 Accredited laboratory). AOCS 0610-A6 is available in 0.5 ml skirted screwcap self-sealing tubes. The soybean FG72 DNA was provided by BASF Agricultural Solutions Seed US LLC and was extracted from clean leaves. The leaf DNA extract sample shall be stored dry in a sealed container at +4 °C in the dark.

Acknowledgements

The authors would like to express sincere appreciation and gratitude to several individuals and their companies for support and guidance throughout this project. Thanks go to Ray Shillito and Benoit Maes, BASF Agricultural Solutions Seed US LLC, for offering AOCS the opportunity to manufacture and distribute these products; to Heather Waxdahl, SGS-Midwest for packaging the samples; and to Dan Smith, Bernd Schoel and Joy Bolster, FoodChain ID, for event-specific, qualitative PCR analysis including the provision of information on running the analyses and interpreting the results.

Glossary

AOCS American Oil Chemists' Society

Conventional Crop A related organism/variety, its components and/or products

for which there is experience of establishing safety based

on common use as food

DNA Deoxyribonucleic Acid is the linear, double-helix

macromolecule that makes up the genetic material of most

organisms

Detection Limit Lowest level at which target DNA can be detected in a sample.

EC European Commission

Genome The full set of genes and associated DNA characteristic of an

organism

GMO Genetically modified/engineered organism: an organism in

which the genetic material has been changed through modern biotechnology in a way that does not occur naturally by

multiplication and/or natural recombination.

ISO International Organisation for Standardisation

PCR Polymerase Chain Reaction: technique used to determine

whether a sample of plant tissue contains a particular DNA

sequence. PCR relies on primer sets that bind to a particular

target DNA sequence and a special DNA-copying enzyme

	(DNA polymerase) that exponentially amplifies the target sequence for identification and measurement
Qualitative PCR	PCR methods that determine the presence or absence of a specific target DNA sequence at a particular level of detection
Quantitation Limit	Lowest level at which the amount of target DNA sequence in a sample can be reliably quantitated
Quantitative PCR	PCR methods that estimate the relative amount of target DNA sequence in a mixture of DNA molecules
Trait: FG72	Tolerance to glyphosate and isoxaflutole herbicides

Introduction

Plant genetic modification is an extension of traditional plant breeding. It allows plant breeders to develop crops with specific traits including insect, disease, and herbicide resistance; processing advantages; and nutritional enhancement. An important component for identifying these new traits is a Certified Reference Material created from leaf, seed, or grain containing the new trait as well as a CRM created from the conventionally bred matrix. The European Commission has mandated that from 18 April 2004, a method for detecting a new event derived from transgenic technology and Certified Reference Material must be available before the EC will consider authorizing acceptance of a new crop derived from transgenic technology. Several nations outside Europe also require grain and ingredients to be labeled above a threshold level before accepting a shipment.

To meet the above regulatory requirements for GMO determination, AOCS 0610-A6 was manufactured from soybean according to ISO 17034:2016 and in accordance with EC No 1829/2003. The CRM is available from AOCS.

Materials and Methods

BASF Agricultural Solutions Seed US LLC prepared the bulk material by taking source leaf material from plants which had been tested individually using several quality standards and was grown from seeds harvested from plants that had themselves passed the same criteria. Plants not meeting the quality standards were removed and destroyed. Leaf material was harvested from the plants which met the quality standards and frozen immediately and stored at -70 °C. The genomic DNA was extracted from leaves of one or more plants according to CTAB-based (Doyle JJ and Doyle JL, 1987) protocol. The integrity and concentration of the genomic DNA was determined by ethidium bromidestained 1% agarose gel electrophoresis and compared to lambda molecular weight standards by digital imaging quantification. The concentration measurement was done in triplicate, repeated in three different gels. No indications for physical degradation were apparent and the DNA migrated at positions higher than 40 Kb.

BASF Agricultural Solutions Seed US LLC delivered 2.0 mg of FG72 soybean leaf DNA

extract to AOCS. The 10 (10) working samples of DNA, 10 µg each, were prepared from

the composite and sent to FoodChain ID, Fairfield, IA (an ISO 17025 Accredited

laboratory) for event-specific, qualitative PCR analysis to screen for the presence of the

intended event, FG72. This testing was for presence confirmation as well as homogeneity

purposes.

The leaf used to manufacture the FG72 materials was shown to contain the FG72 event

as well as the absence of events A2704-12 and A5547-127 (LOQ < 0.1 %) and 3'nos

sequences (LOQ < 0.05 %) using PCR protocols at BASF Agricultural Solutions Seed US

LLC The FG72 soybean leaf DNA was packaged by SGS-Midwest Seed Services in

sterile, 0.5 ml skirted screw-cap self-sealing tubes in aliquots of 10 µg DNA.

AOCS used the Random Number Generator function of Microsoft Excel to select samples

for verification of gene presence, homogeneity, and to rule out degradation during

packaging. Sample numbers AOCS 0610-A6: 14, 44, 52, 61, 76, 95, 115, 141, 155 and

187 were sent to FoodChain ID, Fairfield, IA (an ISO 17025 Accredited laboratory) for

event-specific, qualitative PCR analysis to screen for FG72 presence in the samples.

Stability

Stability of these CRMs has been listed as 1 year from the certification date. The

materials were sealed and stored in the dark at +4 °C, therefore not exposed to air and

are expected to be stable for longer than the estimated expiration date. The stability of

the leaf DNA extract material will be reevaluated annually. If the samples still test positive

for the presence of the trait, the certificates will be extended.

Report of Certification for 0610-A6 Page 9 of 12 ©AOCS, 2022

Results and Discussion

Sample Homogeneity and Prepared Sample Verification

After the bulk material was packaged, ten (10) samples were identified by the Microsoft Excel Random Number Generator and sent to FoodChain ID, Fairfield, IA (an ISO 17025 Accredited laboratory) for event-specific, qualitative PCR analysis. These results are presented in Table 1. This data confirms the presence of the FG72 gene after the packaging of AOCS 0610-A6.

Table 1. Results for the verification of AOCS 0610-A6 FG72 soybean material as tested by FoodChain ID, Fairfield, IA (an ISO 17025 Accredited laboratory) with FG72 event-specific, qualitative PCR analysis.

Sample	FG72 Presence
AOCS 0610-A6 14	Positive
AOCS 0610-A6 44	Positive
AOCS 0610-A6 52	Positive
AOCS 0610-A6 61	Positive
AOCS 0610-A6 76	Positive
AOCS 0610-A6 95	Positive
AOCS 0610-A6 115	Positive

AOCS 0610-A6 141	Positive
AOCS 0610-A6 155	Positive
AOCS 0610-A6 187	Positive

References

Center for Environmental Risk Assessment GM Database http://www.cera-gmc.org/?action=gm_crop_database

FoodChain ID Testing, 504 N. 4th St., Suite 102, Fairfield, IA 52556 Telephone: 1 888 229 2011 www.foodchainid.com

International Seed Testing Association, International Rules of Seed Testing: Seed Science and Technology Rules, 2012

ISO 17025:2005 and ISO 17025:2017, General Requirements for the Competence of Testing and Calibration Laboratories

ISO 174034:2016 (E) General Requirements for the Competence of Reference Material Producers

Regulation (EC) No 1829/2003 of the European Parliament and of the Council of 22 September 2003 on genetically modified food and feed; https://eur-lex.europa.eu/legal-content/EN/TXT/PDF/?uri=CELEX%3A32003R1829&from=en

SGS-Midwest, 236 32nd Avenue, Brookings, SD 57006; Telephone: + 1 605 692 7611 Toll Free: +1 877 692 7611 Fax: +1 605 692 7617 http://www.mwseed.com/