



**Street Address:**  
AOCS, 2710 S. Boulder Drive  
Urbana, IL 61802-6996 USA  
**Phone:** +1-217-359-2344 **Fax:** +1-217-351-8091  
**E-Mail:** CRM@aocs.org **Web:** www.aocs.org

## **Certified Reference Materials**

### **AOCS 0306-G8**

Report of the certification process for

Rf3

Canola Certified Reference Materials

Eighth Batch

Denise Williams  
Technical Services Manager

Scott Bloomer  
Technical Director



ISO 17034:2016  
A2LA Certificate 3438.01

## **Legal Notice**

Neither AOCS nor any person acting on behalf of AOCS is responsible for the use which might be made of the following information.

### **AOCS Mission Statement:**

AOCS advances the science and technology of oils, fats, proteins, surfactants, and related materials, enriching the lives of people everywhere.

More information regarding AOCS is available at <http://www.aocs.org>

## **Contents**

<b>Abstract.....</b>	<b>4</b>
<b>Acknowledgements.....</b>	<b>5</b>
<b>Glossary.....</b>	<b>6</b>
<b>Introduction .....</b>	<b>8</b>
<b>Materials and Methods.....</b>	<b>8</b>
<b>Stability .....</b>	<b>9</b>
<b>Results and Discussion.....</b>	<b>10</b>
<b>Sample Homogeneity and Prepared Sample Verification .....</b>	<b>10</b>
<b>References.....</b>	<b>12</b>

## **Abstract**

This report describes the preparation and certification of the canola CRM AOCS 0306-G8 produced by AOCS Technical Services in 2022. The CRMs have been prepared according to ISO 17034:2016 and are intended to serve as control material for third party testing of canola for transformation events. The presence of Rf3 in the canola was verified using event-specific, qualitative PCR analysis by FoodChain ID Testing, LLC, Chantilly, VA (an ISO 17025 Accredited laboratory). AOCS 0306-G8 is available in 0.5 mL skirted screw-cap self-sealing tubes. The canola Rf3 (breeding line: 0006) DNA was provided by BASF Agricultural Solutions Seed US LLC. The canola Rf3 leaf DNA extract was extracted from clean leaves provided by BASF Agricultural Solutions Seed US LLC. The leaf DNA extract sample shall be stored dry in a sealed container at +4 °C in the dark.

## **Acknowledgements**

The authors would like to express sincere appreciation and gratitude to several individuals and their companies for support and guidance throughout this project. Thanks go to Ray Shillito and Benoit Maes, BASF Agricultural Solutions Seed US LLC, for offering AOCS the opportunity to manufacture and distribute these products; to Heather Waxdahl, SGS-Midwest for packaging the samples; and to Dan Smith, Bernd Schoel and Joy Bolster, FoodChain ID Testing, LLC, for event-specific, qualitative PCR analysis including the provision of information on running the analyses and interpreting the results.

## Glossary

AOCS	American Oil Chemists' Society
Conventional Crop	A related organism/variety, its components and/or products for which there is experience of establishing safety based on common use as food
DNA	Deoxyribonucleic Acid is the linear, double-helix macromolecule that makes up the genetic material of most organisms
Detection Limit	Lowest level at which target DNA can be detected in a sample.
EC	European Commission
Genome	The full set of genes and associated DNA characteristic of an organism
GMO	Genetically modified/engineered organism: an organism in which the genetic material has been changed through modern biotechnology in a way that does not occur naturally by multiplication and/or natural recombination.
ISO	International Organisation for Standardisation
PCR	Polymerase Chain Reaction: technique used to determine whether a sample of plant tissue contains a particular DNA sequence. PCR relies on primer sets that bind to a particular target DNA sequence and a special DNA-copying enzyme

(DNA polymerase) that exponentially amplifies the target sequence for identification and measurement

Qualitative PCR PCR methods that determine the presence or absence of a specific target DNA sequence at a particular level of detection

Quantitation Limit Lowest level at which the amount of target DNA sequence in a sample can be reliably quantitated

Quantitative PCR PCR methods that estimate the relative amount of target DNA sequence in a mixture of DNA molecules

Trait: Rf3 RF3 was produced by genetically engineering plants to restore fertility in the hybrid line and to be tolerant to the herbicide glufosinate ammonium (as a selectable marker).

## **Introduction**

Plant genetic modification is an extension of traditional plant breeding. It allows plant breeders to develop crops with specific traits including insect, disease, and herbicide resistance; processing advantages; and nutritional enhancement. An important component for identifying these new traits is a Certified Reference Material created from leaf, seed, or grain containing the new trait as well as a CRM created from the conventionally bred matrix. The European Commission has mandated that from 18 April 2004, a method for detecting a new event derived from transgenic technology and Certified Reference Material must be available before the EC will consider authorizing acceptance of a new crop derived from transgenic technology. Several nations outside Europe also require grain and ingredients to be labeled above a threshold level before accepting a shipment.

To meet the above regulatory requirements for GMO determination, AOCS 0306-G8 was manufactured from canola according to ISO 17034:2016 and in accordance with EC No 1829/2003. The CRM is available from AOCS.

## **Materials and Methods**

BASF Agricultural Solutions Seed US LLC prepared the bulk material by taking source leaf material from plants which had been tested individually using several quality standards and was grown from seeds harvested from plants that had themselves passed the same criteria. Plants not meeting the quality standards were removed and destroyed. Leaf material was harvested from the plants which met the quality standards and frozen immediately and stored at -70°C. The genomic DNA was extracted from leaves of one or more plants according to CTAB-based (Doyle JJ and Doyle JL, 1987) protocol. The integrity and concentration of the genomic DNA was determined by ethidium bromide-stained 1% agarose gel electrophoresis and compared to lambda molecular weight standards by digital imaging quantification. The concentration measurement was done in triplicate, repeated in three different gels. No indications for physical degradation were apparent and the DNA migrated at positions higher than 40 Kb.



BASF Agricultural Solutions Seed US LLC delivered 3.0 mg of Rf3 canola leaf DNA extract to AOCS. The ten (10) working samples of DNA, 10 µg each, were prepared from the composite and sent to FoodChain ID Testing, LLC, Chantilly, VA (an ISO 17025 Accredited laboratory) for event-specific, qualitative PCR analysis to screen for the presence of the intended event, Rf3. This testing was for presence confirmation as well as homogeneity purposes.

The leaf used to manufacture the Rf3 materials was shown to contain the Rf3 event as well as the absence of 2mEPSPS, cp4EPSPS, BXN, nptII, Barnase, pat and P35S (LOD < 0.05%) using PCR protocols at BASF Agricultural Solutions Seed US LLC. The Rf3 canola leaf DNA was packaged by SGS-Midwest Seed Services in sterile, 0.5 mL skirted screw-cap self-sealing tubes in aliquots of 10 µg DNA.

AOCS used the Random Number Generator function of Microsoft Excel to select samples for verification of gene presence, homogeneity, and to rule out degradation during packaging. Sample numbers AOCS 0306-G8: 13, 38, 61, 70, 82, 89, 100, 125, 265 and 295 were sent to FoodChain ID Testing, LLC, Chantilly, VA (an ISO 17025 Accredited laboratory) for event-specific, qualitative PCR analysis to screen for Rf3 presence in the samples.

## **Stability**

Stability of these CRMs has been listed as 1 year from the certification date. The materials were sealed and stored in the dark at 4 °C, therefore not exposed to air and are expected to be stable for longer than the estimated expiration date. The stability of the leaf DNA extract material will be reevaluated annually. If the samples still test positive for the presence of the trait, the certificates will be extended.

## Results and Discussion

### Sample Homogeneity and Prepared Sample Verification

After the bulk material was packaged, ten (10) samples were identified by the Microsoft Excel Random Number Generator and sent to FoodChain ID Testing, LLC, Chantilly, VA (an ISO 17025 Accredited laboratory) for event-specific, qualitative PCR analysis. These results are presented in Table 1. This data confirms the presence and homogeneity of the Rf3 gene after the packaging of AOCS 0306-G8.

Sample	Rf3 Presence
AOCS 0306-G8 13	Positive
AOCS 0306-G8 38	Positive
AOCS 0306-G8 61	Positive
AOCS 0306-G8 70	Positive
AOCS 0306-G8 82	Positive
AOCS 0306-G8 89	Positive
AOCS 0306-G8 100	Positive

AOCS 0306-G8 112	Positive
AOCS 0306-G8 265	Positive
AOCS 0306-G8 295	Positive

## References

Center for Environmental Risk Assessment GM Database

[http://www.cera-gmc.org/?action=gm\\_crop\\_database](http://www.cera-gmc.org/?action=gm_crop_database)

FoodChain ID Testing, LLC, 4150 Lafayette Center Drive, Suite 600, Chantilly, VA 20151  
www.foodchainid.com

International Seed Testing Association, International Rules of Seed Testing: Seed Science and Technology Rules, 2012

ISO 17025:2005 and ISO 17025:2017, General Requirements for the Competence of Testing and Calibration Laboratories

ISO 174034:2016 (E) General Requirements for the Competence of Reference Material Producers

Regulation (EC) No 1829/2003 of the European Parliament and of the Council of 22 September 2003 on genetically modified food and feed; <https://eur-lex.europa.eu/legal-content/EN/TXT/PDF/?uri=CELEX%3A32003R1829&from=en>

SGS-Midwest, 236 32nd Avenue, Brookings, South Dakota 57006 Telephone: + 1 605 692 7611 Toll Free: +1 877 692 7611 Fax: +1 605 692 7617 <http://www.mwseed.com/>